



LUVATOL

Thickening Pastes for SMC / BMC

Productname / Thickening			MgO	Viscosity	Density	Application				
medium thickening	high thickening	low thickening	% Active Ingredient	mPas (25° C)	g/ccm (25° C)	Standard	Low Shrink	Low Profile	Low Pressure	Relining
Luvatol MK 15 N			15	4.000	1.50	X	X	X	X	
Luvatol MK 25			25	4.000	1.55	X	X	X	X	X
Luvatol MK 2050			25	20.000	1.75	X	X	X	X	
Luvatol MK 25 MZ			25	10.000	1.68	X	X	X	X	
Luvatol MK 25 NV			25	1.000	1.40	X	X	X	X	X
Luvatol MK 35			35	17.000	1.71	X	X	X		
Luvatol MK 35 NV			35	10.000	1.68	X	X	X		
Luvatol MK 36			35	11.000	1.67	X	X	X		
Luvatol MK 40			40	10.000	1.65	X	X	X		
Luvatol EK 25 NV			25	4.500	1.55	X	X	X	X	X
Luvatol EK 30 NV			30	200	1.21	X	x	X	X	X
Luvatol EK 35 KM			35	6.000	1.57	X	X	X		
Luvatol EK 100 KM			50	8.000	1.60	X	X	X		
Luvatol DH 25			27	8.000	1.55			X	X	
Luvatol EH 35			35	11.000	1.58			X	X	
Luvatol DK 35			35	8.000	1.55			X	X	
Luvatol DC 35			35	4.000	1.55					
Luvatol HP 40			40	8.000	1.55		X	X	X	

Issue: 02/2014

Contact: luvaton@lehvoss.de